



# STERIPHENE II®

brand disinfectant deodorant

Steriphene II Brand Disinfectant Deodorant is a versatile, ready to use aerosol disinfectant deodorant. It is ideal for convenient disinfection after spot clean up, and for those areas not easily reached during routine disinfection procedures. Because it's an aerosol, there are no mixing problems with Steriphene II nor any waste.



Providing Quality Maintenance and Chemical Specialty Solutions with Unparalleled Customer Service



Spartan Chemical Company, Inc.  
 1110 Spartan Drive  
 Maumee, OH 43537  
 1-800-537-8990  
 www.spartanchemical.com

**HOSPITAL TYPE DISINFECTANT:**

Steriphene II is a “hospital type disinfectant” which achieves 100% kill of all claimed microorganisms (see product label). A blend of two phenolics, Steriphene II is Tuberculocidal, Bactericidal, Fungicidal\* and Virucidal\*\* and is effective against Human Immunodeficiency Virus Type 1 (AIDS Virus), HBV (Hepatitis B Virus), Herpes Simplex Type 2 (causes genital herpes) and Influenza A/J305. EPA Reg. No. 5741-22.

**EFFECTIVE AGAINST ANTIBIOTIC RESISTANT BACTERIA\*\*\*:**

Steriphene II has been tested and shown effective against Methicillin-resistant Staphylococcus aureus (MRSA) and Vancomycin-resistant Enterococcus faecalis (VRE).

**EXCELLENT DEODORIZER:**

Steriphene II destroys bacteria the initial source of odors. It is not a cover up. Deodorizers which do not have disinfectant properties merely cover one odor with another...in effect, desensitizing the sense of smell.

**CHOOSE ONE OF TWO UNIQUE FRAGRANCES!**

Steriphene II is offered in two distinctive fragrances...Spring Breeze and Clean, Fresh.

**HIGHLY ACTIVE:**

Steriphene II’s active ingredient content is 64.122%...as much as 1/3 more active ingredients than many competitive products.

**FORMERLY AUTHORIZED BY THE USDA:**

Category “C1” – compound for use on all surfaces in inedible product areas, non-processing areas, and/or exterior areas.

**DIRECTIONS FOR USE:**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Clean all surfaces prior to disinfection.

**Disinfection and Deodorization:** Shake well. Hold can 6 to 8 in. from surface and apply in a uniform moist film. Treated surfaces should remain wet for 10 min. Allow to air dry or wipe off.

**Tuberculocidal Activity:** Steriphene II Brand Disinfectant is tuberculocidal under the above conditions in 10 min. at 68°F (20°C).

**\*Kills HBV and HIV on pre-cleaned environmental surfaces/ objects previously soiled with blood/body fluids** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Hepatitis B virus (HBV) and human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

**Special instructions for cleaning and decontamination against HBV and HIV-1 on surfaces or objects soiled with blood/body fluids:**

- 1. Personal Protection:** Clean up should be done wearing protective latex gloves, gowns, masks and eye protection.
- 2. Cleaning Procedure:** Blood or body fluids must be thoroughly cleaned from surfaces and objects before application of disinfectant.

- 3. Contact Time:** Leave surfaces wet for 10 minutes.
- 4. Disposal of Infectious Material:** Contaminated cleaning materials, blood and other body fluids should be autoclaved and/or disposed of according to local regulations for infectious waste disposal.

**NOTE:** Make sure valve opening is pointed in the direction of intended spray. Do not use on varnished wood or rayon fabrics. Steriphene II Brand Disinfectant Deodorant deodorizes at the source by killing the germs that cause malodors. This product is not intended for use as a general space spray air freshener.

**Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.**

**SPECIFICATION DATA:**

Contains sodium nitrite  
 pH - 8.0-9.0  
 Fragrances - Clean Fresh or Spring Breeze  
 Flash Point - Flammable  
 Stability - Shelf @ 24°C/75°F - One year minimum  
 Contains biodegradable surfactants  
 Former USDA rating: C1  
 Does not require Kosher certification per Organized Kashruth Laboratories.



Distributed by:



1890 Chrysler Dr.  
 Belvidere, IL 61008  
 815-547-3737  
 sales@2hdi.com  
 www.HDi3.com

PACKAGING		
608100	20-ounce cans (12 per case)	Clean Fresh
607500	20-ounce cans (12 per case)	Spring Breeze



Label copy is provided in English.

GUARANTEE: Spartan’s modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 5/11 L6081



# SparSan Q<sup>®</sup>

## disinfectant deodorant

SparSan Q Disinfectant Deodorant is a versatile, ready to use aerosol disinfectant deodorant. It is ideal for convenient disinfection after spot clean up, and for those areas not easily reached during routine disinfection procedures. Because it's an aerosol, there is no mixing or any waste.





Spartan Chemical Company, Inc.  
 1110 Spartan Drive  
 Maumee, OH 43537  
 1-800-537-8990  
 www.spartanchemical.com

- Kills HIV-1 on Pre-Cleaned Environmental Surfaces
- Eliminates odors and has a linen clean fragrance
- Germicidal
- Fungicidal
- Pseudomonacidal
- Tuberculocidal
- Kills Household Germs
- Prevents Mold and Mildew
- Kills Methicillin-Resistant Staphylococcus Aureus (MRSA), Vancomycin-Resistant Enterococcus Faecalis (VRE), Influenza A Virus, Influenza B Virus, Swine Influenza A Virus Strain (H1N1), and Influenza A Virus Strain Hong Kong
- **Sanitizes:** Kills 99.9% of Staphylococcus aureus (Staph), and Enterobacter aerogenes on pre-cleaned, hard nonporous, non-food contact surfaces with a 5 minute contact time.

#### HOSPITAL-GRADE DISINFECTANT

Disinfects and deodorizes surfaces through effective bactericidal and fungicidal activity. Kills Staphylococcus aureus (Staph), Streptococcus hemolyticus (Strep), Pseudomonas aeruginosa, Mycobacterium bovis (BCG), Salmonella choleraesuis and Trichophyton mentagrophytes (athletes foot fungus) on hard, non-porous environmental surfaces. Meets EPA Efficacy Standards for hospital aerosol disinfectants.

#### SPRAYED ON SURFACES

Deodorizes by disinfecting. Germicidal activity kills most odor-producing organisms. Controls unpleasant odors. Pleasantly scented. For use in hospital rooms, sick rooms, doctors' examining rooms, kitchens, nursery and bathrooms.

#### EXCELLENT DEODORIZER

SparSan Q Disinfectant Deodorant destroys bacteria the initial source of odors. It is not a cover up. Deodorizers which do not have disinfectant properties merely cover one odor with another.



Distributed by:



1890 Chrysler Dr.  
 Belvidere, IL 61008  
 815-547-3737  
 sales@2hdi.com  
 www.HDi3.com

#### DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Shake well before using. Prior to application, thoroughly clean all surfaces to be treated to remove gross filth or heavy soil.

**To Disinfect:** Hold can upright and spray from a distance of 6 to 8 inches until area for treatment is thoroughly moistened and allow to dry for 10 minutes before wiping or allow to air dry.

**For Control of Mold and Mildew on Hard, Non-Porous Environmental Surfaces:** Spray surface until thoroughly moistened. Allow to air dry. Repeat application in 7 days.

**For Control of Fungus:** In areas such as bathrooms, locker rooms, floors, benches, shower stalls or other non-porous environmental surfaces commonly contacted by bare feet: Spray surface until thoroughly moistened. Allow to air dry.

**For Use On:** Walls, floors, bed frames, bedpans, urinals, fixtures, containers, toilet seats, diaper pails, telephones, trash cans, and sickrooms. Due to the possibility of electric shock when used on telephones, avoid excess spraying of telephone.

*Note: Disinfected food contact surfaces must be rinsed with potable water after cleaning and/or disinfection and before reuse.*

#### Kills Hiv On Pre-Cleaned Environmental Surfaces And Objects Previously Soiled With Blood/Body Fluids:

In health care settings. Or in other settings where there is a likelihood of soiling of inanimate surfaces or objects with blood or body fluids in which the surfaces and objects can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

#### Special instructions for cleaning and decontamination against hiv on surfaces and objects soiled with blood or body fluids:

##### Personal protection:

Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks or eye coverings.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of disinfectant.

**Disposal Of Infectious Material:** Blood and other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

**Contact Time:** Leave surfaces wet for 10 minutes. This contact time will not control other common types of viruses and bacteria.

**Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting [www.spartanchemical.com](http://www.spartanchemical.com).**

#### ACTIVE INGREDIENTS

n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	
dimethyl benzyl ammonium chlorides...	00.05%
n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides.....	00.05%
Ethyl Alcohol.....	37.00%
<b>OTHER INGREDIENTS.....</b>	<b>62.90%</b>
<b>Total.....</b>	<b>100.00%</b>

*Contains Sodium Nitrite*

#### SPECIFICATION DATA

**Form:** Pressurized liquid (aerosol can)  
**pH:** 10.5-11.5  
**Flash Point:** Flammable  
**Stability:** Shelf @ 24°C/75°F: One year minimum



#### PACKAGING

**607600** 20 Ounce Cans (12 Per Case)

**GUARANTEE:** Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 09/15 L6076



# HARD SURFACE DISINFECTING WIPES

## NABC®, Fresh, and Lemon scents



Cross-contamination of hard surfaces happens anywhere people work and play. Assisted living facilities and nursing homes, doctors' offices, daycare centers, grocery stores, and fitness centers are just a few of the locations where viruses, bacteria, and disease can be spread from one hard surface to another. Reducing the spread of germs must be part of a daily routine. The solution: wipe it away!



1890 Chrysler Dr.  
Belvidere, IL 61008  
815-547-3737  
sales@2hdi.com  
www.HDi3.com

*Providing quality maintenance and chemical specialty solutions with unparalleled customer service since 1956*



# HARD SURFACE DISINFECTING WIPES

## FEATURES

- Cleans, disinfects<sup>3</sup>, and deodorizes
- Kills 99.9% of bacteria in just 15 seconds
- Kills the cold<sup>1</sup> & flu viruses on hard non-porous surfaces<sup>†††</sup>
- Available in Lemon, Fresh, and NABC scents
- Use in healthcare settings, schools, grocery, athletic, and other public facilities

## BENEFITS

- Convenient and ready to use for cleaning and disinfecting hard, non-porous surfaces
- Fast and effective cleaning
- Effective against odor-causing bacteria

<sup>1</sup> This product is an effective virucide on hard non-porous surfaces against Human Coronavirus (ATCC VR-740) when the treated surface is allowed to remain wet for 10 minutes. Allow surfaces to air dry.

<sup>3</sup> When used according to disinfection directions

<sup>†††</sup> Influenza A Virus, Pandemic 2009 H1N1 Influenza A Virus, Human Coronavirus<sup>1</sup>

## SPECIFICATION DATA

pH	10.5 – 12.0
Fragrance	NABC <sup>®</sup> Scent Fresh Scent Lemon Scent
Specific Gravity @ 24°C/75°F	1.007
EPA Registration Number	1839-190-5741

## RELATED PRODUCTS



FLOOR  
STAND  
#930700



STAND  
EXTENSION  
#930701



WALL  
BRACKET  
#930800

\* Each canister contains 125 pre-moistened, 7" x 8" wipes.

**SAFETY:** Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at [www.spartanchemical.com](http://www.spartanchemical.com).

**GUARANTEE:** Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.

## EFFICACY DATA

### DISINFECTING

This product is an effective disinfectant against *Escherichia coli* O157:H7, *Pseudomonas aeruginosa*, *Salmonella enterica*, *Staphylococcus aureus*, *Streptococcus pyogenes*, Vancomycin resistant *Enterococcus faecium*, Methicillin resistant *Staphylococcus aureus*, Community Associated Methicillin resistant *Staphylococcus aureus* Genotype USA400, Community Associated Methicillin resistant *Staphylococcus aureus* Genotype USA300 and Vancomycin intermediate resistant *Staphylococcus aureus* and *Klebsiella pneumoniae* New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant (CDC Clinical Isolate No. 10002).

### VIRUCIDAL ACTIVITY

This product is an effective virucide on hard, non-porous surfaces against Influenza A2/Japan (VR-547), Herpes Simplex Virus 1 (ATCC VR- 260), Human Immunodeficiency Virus, Type 1 (HIV-1) and Herpes Simplex Virus 2 (ATCC VR-734), Duck Hepatitis B Virus (DHBV), Human Hepatitis B Virus (HBV), Bovine Viral Diarrhea Virus (BVDV) and Human Hepatitis C Virus (HCV) when the treated surface is allowed to remain wet for 5 minutes. Allow surface to air dry.

<sup>1</sup> This product is an effective virucide on hard non-porous surfaces against Human Coronavirus (ATCC VR-740) when the treated surface is allowed to remain wet for 10 minutes. Allow surfaces to air dry.

**KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1, Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

### FUNGICIDAL ACTIVITY

This product is an effective fungicide against *Trichophyton interdigitale*, formerly *Trichophyton mentagrophytes* (causative agent of Athlete's Foot fungus) (a cause of ringworm) when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 3 minutes. Allow surfaces to air dry.

### SANITIZING/DEODORIZING

To sanitize/deodorize hard, non-porous surfaces: wipe surface, using enough wipes to allow treated surfaces to remain wet for 15 seconds for *Staphylococcus aureus*, *Klebsiella pneumoniae*, *Listeria monocytogenes*, *Escherichia coli* O157:H7 and 60 seconds for *Campylobacter jejuni*. Let surface dry. For highly soiled surfaces, clean excess dirt first.

## ORDERING INFORMATION\*



#108506  
Lemon  
Scent



#108606  
Fresh  
Scent



#108706  
NABC<sup>®</sup>  
Scent





# PROFECT® HP

## Hydrogen Peroxide Disinfectant



Use Profect HP daily as part of a simple and effective cleaning and disinfection program for your entire facility on hard, non-porous surfaces. Featuring patented hydrogen peroxide technology, Profect HP kills bacteria and viruses in 60 seconds. One minute contact times ensure efficacy and compliance for your most critical disinfection needs. Ideal for daily use on high-touch surfaces, Profect HP is available in a convenient ready-to-use quart.

EPA Registration Number: 5741-33

# PROFECT® HP

## FEATURES

- Kills bacteria, enveloped and non-enveloped viruses, and fungi
- Features patented hydrogen peroxide technology
- Ready-to-use formula ensures proper ppm required for disinfection
- Kills SARS-CoV-2 in 30 seconds

## BENEFITS


- One minute contact times ensure efficacy and compliance
- Offers both safety and broad-spectrum kill in one product
- Meets the surface disinfection recommendations of OSHA bloodborne pathogen standard
- Disinfects high-touch surfaces in 1 minute

## SPECIFICATION DATA

Dilution	Ready-to-use
pH	1.8 – 2.2
Color	Colorless
Fragrance	Mild
Specific Gravity @ 24°C/75°F	1.014
Flash Point	> 94°C/201°F ASTM D56
Stability: Shelf @ 24°C/75°F	1 year minimum
Stability: Freeze/Thaw	Will withstand 3 cycles
EPA Registration Number	5741-33

## ADDITIONAL RESOURCES

Scan the QR code to access this product's safety information and documentation.



## EFFICACY DATA

### FOR USE AS A ONE-STEP DISINFECTANT, DEODORIZER, AND CLEANER ON HARD, INANIMATE NONPOROUS ENVIRONMENTAL SURFACES

When used as directed, this product is effective against *Pseudomonas aeruginosa*, *Salmonella enterica*, *Staphylococcus aureus*, *Acinetobacter baumannii*, *Escherichia coli* 0157:H7, Methicillin Resistant *Staphylococcus aureus* (MRSA), Multi-drug Resistant (MDR) *Enterococcus faecium*.

### TUBERCULOCIDAL ACTIVITY

Profect HP exhibits disinfectant efficacy against *Mycobacterium bovis* (TB) when the treated surface is allowed to remain wet for 3 minutes at 68°F (20°C).

### FUNGICIDAL ACTIVITY

Profect HP is an effective disinfectant against the pathogenic fungus, *Trichophyton interdigitale* (Athlete's Foot Fungus) (a cause of ringworm) on visibly clean, hard, non-porous surfaces (such as those found in locker rooms, dressing rooms, shower and bath areas, and exercise areas) for a 1-minute contact time at room temperature.

### \*VIRUCIDAL ACTIVITY ON HARD, INANIMATE NONPOROUS ENVIRONMENTAL SURFACES

When used as directed, Profect HP is effective against Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), Human Immunodeficiency Virus Type 1 (HIV-1), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Coronavirus, Human Rotavirus, Influenza A Virus, Influenza B Virus, and Norovirus when treated surface is allowed to remain wet for 1 minute.

When used as directed, Profect HP is effective against Severe Acute Respiratory Syndrome-related Coronavirus 2 (SARS-COV-2) in 30 seconds!

## ORDERING INFORMATION



#100803



1890 Chrysler Dr.  
Belvidere, IL 61008  
815-547-3737  
sales@2hdi.com  
www.HDi3.com



**SAFETY:** For institutional and industrial use only. Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at [www.spartanchemical.com](http://www.spartanchemical.com). Workplace labels are available in English, Spanish, and French.

**GUARANTEE:** Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.

[www.spartanchemical.com](http://www.spartanchemical.com) • 1-800-537-8990 • 1110 Spartan Drive, Maumee, OH 43537

© 2021 Spartan Chemical Company, Inc. All Rights Reserved. Profect is a registered trademark of Spartan Chemical Company, Inc.



11/21 L1008