# INFORMATIVE BULLETIN



Spartan Chemical Company, Inc., 1110 Spartan Drive, Maumee, OH 43537 Toll Free: 1-800-537-8990, Fax: 419-536-8423, Fax-On-Demand: 888-383-4666 E-Mail: spartan@spartanchemical.com, Web Site: www.spartanchemical.com @SCC 2003







The Carpet and Rug Institute's Seal of Approval signifies this product effectively removes even the most difficult stains without any damage to the carpet.

#### PRODUCT DESCRIPTION:

Contempo  $H_2O_2$  Spotting Solution is a convenient, ready-to-use, hydrogen peroxide- based carpet spotting solution and pre-spray treatment. The clear colorless solution has a fresh crystal water fragrance.

### **EFFERVESCENT BUBBLES = CLEANING ACTION:**

Contempo  $H_2O_2$  Spotting Solution provides *on-the-spot* cleaning action. This hydrogen peroxide-based product cleans away many of the most dreaded food stains: red wine, red sauces, mustard, grape juice, coffee, tea, food grade red dye, and chocolate. Contempo  $H_2O_2$  is also effective against most cooking oils and greases, including motor oil. The incorporation of hydrogen peroxide in Contempo  $H_2O_2$  Spotting Solution provides the unique ability to oxidize odors from stains caused by pet soils, vomit, and spoiled milk.

#### NO TACKY RESIDUE:

Contempo  $H_2O_2$  Spotting Solution is an advanced cleaning formula designed with ingredients that penetrate, lift and oxidize stains. These powerful detergents dry to a brittle powder, leaving no tacky residue. The brittle powder will not contribute to re-soiling and is easily removed by vacuuming.

### SPOT VS. STAIN PHILOSOPHY:

Spots that are left untreated are future stains. Quick, timely action with this superior spotting agent is the best defense for spot removal and stain prevention.

#### **DIRECTIONS FOR USE:**

Some individuals may be sensitive to ingredients in this product. Before using, read product label and MSD sheet. Wet floors may be slippery. Prevent pedestrian traffic with signs or barricades.

Before use, always pre-test product on an inconspicuous area to test for colorfastness, bleaching or shrinkage.



#### Spotting:

- 1. Remove excess matter and blot up excess moisture with a clean, white, absorbent cloth.
- 2. Spray Contempo  $H_2O_2$  Spotting Solution directly on the spot and wait a few seconds to allow the product time to penetrate.
- 3. Gentle agitation may be needed. Use a carpet taping knife or soft carpet brush to penetrate the carpet pile, working the cleaning solution into the affected area.
- 4. Spray a small amount of water on the area to float the spot up and out of the carpet pile.
- 5. Blot with a clean, white, absorbent towel. Frequently turn the towel and apply a clean area of the cloth to the carpet until the spot is lifted. Extraction is recommended if available.
- 6. If the stain is set deeply into the pile, repeated applications may be necessary. Since liquids will spread as they penetrate the nap, you will need to spray a surface area twice as large as the stain. A wet-vac pick-up will help speed the drying process. Vigilant attention to spots can extend the time between all-over carpet cleanings.

## **Pre-Spray Extract Method:**

- 1. Apply Contempo H<sub>2</sub>O<sub>2</sub> Spotting Solution directly from the product container or through a pump-up or pressure type applicator prior to shampooing or extracting. Do not dilute. (Spartan Defoamer may be needed in the recovery tank.)
- 2. Let set 15 minutes.
- 3. Extract with clear water.
- 4. When carpet is completely dry, vacuum thoroughly.

## **Pre-Treating Traffic Lanes and Other Heavily Soiled Areas:**

Apply directly to carpet from the product container or through a pump-up or pressure type applicator prior to shampooing or extracting. Do not dilute. Allow product to sit and penetrate the carpet pile. Shampoo or extract as usual.

#### **SPECIFICATION DATA:**

Specific Gravity – 1.016 @ 24°C/75°F pH – 5.0 – 6.0 Density – 8.46 lbs/gal @ 24°C/75°F Shelf

- a. Shelf @ 24°C/75°F One year
- b. Freeze/Thaw Cycle 3 cycles minimum

#### **PACKAGING:**

Contempo H<sub>2</sub>O<sub>2</sub> Spotting Solution is available in carafe-style quart bottles, 12 per case, and gallons, 4 per case. Each case of quart containers includes a pair of gloves and three trigger sprayers. Each case of gallons includes a pair of gloves. Label copy is provided in English, Spanish, and French.

Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.

#### **GUARANTEE:**

Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture.

